



# Recombinant Danio rerio Paired box protein Pax-6 (pax6a)

<b>Product Code</b>	CSB-MP017492DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P26630
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPQKEYYNRA TWESGVASMM QNSHSGVNQL GGVFVNGRPL PDSTRQKIVE LAHSGARPCD ISRILQVSNG CVSKILGRYY ETGSIRPRAI GGSKPRVATP EVVGKIAQYK RECPSIFAW E IRDRLLSEGV CTNDNIPSVS SINRVLRLNLA SEKQQMGADG MYEKLRLMNG QTGTWGTTPG WYPGTSVPGQ PNQDGCQQSD GGGENTNSIS SNGEDSDETQ MRLQLKRKLQ RNRTSFTQEQ IEALEKEFER THYPDVFARE RLA AKIDLPE ARIQVWFSNR RAKWRREEKL RNQRRQASNS SSHIPISSSF STSVYQPIPQ PTTPVSFTSG SMLGRSDTAL TNTYSALPPM PSFTMANNLP MQPSQTSSYS CMLPTSPSVN GRSYDITYTPP HMQAHMNSQS MAASGTTSTG LISPGVSVPV QVPGSEPDMS QYWPRLQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pax6a
<b>Protein Names</b>	Recommended name: Paired box protein Pax-6 Alternative name(s): Pax[Zf-a]
<b>Expression Region</b>	1-437
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.