



Recombinant Sheep Plasminogen (PLG)

Product Code	CSB-BP018188SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81286
Product Type	Recombinant Protein
Immunogen Species	Ovis aries (Sheep)
Purity	>85% (SDS-PAGE)
Sequence	APQAPSVENP PEADCMLGIG KGYRGKKATT VAGVPCQEWA AQEPHRHGIF TPETNPRAGL EKNYCRNPDG DVNGPWCYTT NPRKLFDYCD IPQCESSFDC GKPKVEPKKC PARVVGGCVA TPHSWPWQVS LRRRSREHFC GGTLISPEWV LTAAHCLDSI LGPSFYTVIL GAHYEMAREA SVQEIPVSRL FLEPSRADIA LKLSPPAVI TDEVIPACLP SPNYVVADKT VCYITGWGET QGTFGVGRLK EARLPVIENK VCNRYEYLNQ RVKSTELCAG DLAGGTDSCQ GDSGGPLVCF EKDKYILQGV TSWGLGCARP NKPGVYVRVS TYVPWIEETM RRY
Source	Baculovirus
Target Names	PLG
Protein Names	Recommended name: Plasminogen EC= 3.4.21.7
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.