



Recombinant *Saccharomyces cerevisiae* Diphthine synthase (DPH5)

Product Code	CSB-MP007133SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32469
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MLYLIGLGLS YKSDITVRGL EAIKKCSRVY LEHYTSILMA ASQEELESYY GKEILADRE LVETGSKQIL NNADKEDVAF LVVGDPFGAT THTDLVLRAK REAIPVEIIH NASVMNAVGA CGLQLYNFGQ TVSMVFFTDN WRPDSWYDKI WENRKIGLHT LVLLDIKVKE QSIENMARGR LIYEPPRYMS IAQCCEQLLE IEEKRGTKAY TPDTPAV AIS RLGSSSQSFK SGTISELANY DSGEPLHSLV ILGRQCHELE LEYLLEFADD KEKFGKDVAN DQEYFKPAAW VPPTEDDSDE
Source	Mammalian cell
Target Names	DPH5
Protein Names	Recommended name: Diphthine synthase EC= 2.1.1.98 Alternative name(s): Diphthamide biosynthesis methyltransferase
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.