



Recombinant Mouse Pentraxin-related protein PTX3 (Ptx3)

| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP019079MO-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P48759 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | ETS DDYELMYVNL DNEIDNGLHP TEDPTPCDCR QEHSEWDKLF IMLENSQMRE GMLLQATDDV LRGELQRLRA ELGRLAGGMA RPCAAGGPAD ARLVRALEPL LQESRDASLR LARLEDAEAR RPEATVPGLG AVLEELRRTR ADLSAVQSWV ARHWLPAGCE TAIFFPMSK KIFGSVHPVR PMKLESFSTC IWKATDVLN KTILFSYGTK WNPYEIQLYL SSQSLVLVVG GKENKLAADT VVSLGRWSHL CGTWSSEQGS MSLWANGELV ATTVEMAKSH SVPEGGLLQI GQEKNGCCVG GGFDESLAFS GRITGFNIWD RVLSEEEIRA SGGVESCHIR GNVVGWGVTE IQAHGGAQYV S |
| Source | E.coli |
| Target Names | Ptx3 |
| Protein Names | Recommended name: Pentraxin-related protein PTX3 Alternative name(s): Pentaxin-related protein PTX3 Tumor necrosis factor-inducible gene 14 protein Short name= TSG-14 |
| Expression Region | 18-381 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |