



Recombinant *Saccharomyces cerevisiae* Porphobilinogen deaminase (HEM3)

Product Code	CSB-EP010524SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28789
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGPETLHIGG RKSKLAVIQS NHVLKLIIEK YPDYDCKVFT LQTLGDQIQF KPLYSGGKA LWTKELEDHL YHDDPSKKLD LIVHSLKDMP TLLPEGFELG GITKRVDPTD CLVMPFYSAY KSLDDLDPGG IVGTSSVRRS AQLKRKYPHL KFESVRGNIQ TRLQKLDDPK SPYQCIILAS AGLMRMGLEN RITQRFHSDT MYHAVGQGAL GIEIRKGDTK MMKILDEICD LNATICCLSE RALMRTLEGG CSVPIGVESK YNEETKLLLL KAIVVDVEGT EAVEDEIEML IENVKEDSMA CGKILAERMI ADGAKKILDE INLDRIK
Source	E.coli
Target Names	HEM3
Protein Names	Recommended name: Porphobilinogen deaminase Short name= PBG EC= 2.5.1.61 Alternative name(s): Hydroxymethylbilane synthase Short name= HMBS Pre-uroporphyrinogen synthase
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.