



Recombinant Dictyostelium discoideum Tubulin beta chain (tubB)

Product Code	CSB-MP025318DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32256
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MREIVQIQAG QCGNQIGSKF WEVISEEHGI QSDGFHAGGE DEHLKRLQLE RINVYYNEAR DGKYVPRSVL VDLEPGTVDI IKASPYGKLF RPDNFIHGQS GAGNNWAKGH YTEGVVELVES VLDVVRRETE GCDCLQGFQV THSIGGGTGS GLGTLISKI REEFPDRMMC TFSVVPSPKV SLTVVEPYNA TLVHQLVEN ADEVMCIDNE ALHDICFRTL KLTQPNYGD L NLISSVMSG ITCCLRFPGQ LNSDLRKLAV NLIPFRLHF FLVGFAPLTA KGASSYNRIT VPELTQQMFD AKNMMAASDP HNGKYL TASA LFRGKIFTKE VDEQMHN IQT KNSSYFVEWI PHNIKSSICD IPPKGT PMAV TFIGNNTAIQ DLFKRIS IHF QAMFRRKAFL HWYTLEG MEE LEFTEAESNM NDLVYEQY QY SNQETEEDGG EYQEEHEEHE EQAEN
Source	Mammalian cell
Target Names	tubB
Protein Names	Recommended name: Tubulin beta chain Alternative name(s): Beta-tubulin
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.