



Recombinant Putative N-acetylglucosamine-6-phosphate deacetylase (F59B2.3)

Product Code	CSB-BP001660CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34480
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLPRCLKIDS TELAEITNEL VQFLNCLVLR SGGLKKEHIW VRNGRILDER TVFFEEKTMA DVQIDCEGLI LSPGFIDLQL NGGFGIDFST YNSDDKEYQE GLALVAKQLL AHGVTFSFPT VITSSPETYH KILPLLKPSN ASSEGAGNLG AHLEGPFISA DKRGCHPEQL VITSLSPNPV EIIEHVYGST ENIAIVTMAP ELEGAQEAIE YFVSTGTTVS VGHSSAKLGP GEMAVLSGAK MITHLFNAMQ SYHHRDPGLI GLLTSSKLTP DHPLYYYGIIS DGIHTHDSAL RIAYHTNSAG LVLVTDIAAA LGMSDGVHKL GTQTIHVKGL EAKLDGTNTT AGSVASMPYC IRHLMKATGC PIEFALQSAT HKPATLLGVS DEKGTLDVGR LADFVLIDKN VTVKATFCSG KRVFLAQD
Source	Baculovirus
Target Names	F59B2.3
Protein Names	Recommended name: Putative N-acetylglucosamine-6-phosphate deacetylase EC= 3.5.1.25 Alternative name(s): GlcNAc 6-P deacetylase
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.