



# Recombinant Human Deoxyribonuclease-1-like 1 (DNASE1L1)

<b>Product Code</b>	CSB-MP007050HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P49184
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	FR ICAFNAQRLT LAKVAREQVM DTLVRILARC DIMVLQEVVD SSGSAIPLLL RELNRFDGSG PYSTLSSPQL GRSTYMETYV YFYRSHKTQV LSSYVYNDED DVFAREFPVA QFSLPSNVLP SLVLVPLHTT PKAVEKELNA LYDVFLEVSQ HWQSKDVILL GDFNADCASL TTKRLDKLEL RTEPGFHWVI ADGEDTTVRA STHCTYDRVV LHGERCRSL LHTAAAFDFPT SFQLTEEEAL NISDHYPVEV ELKLSQAHSV QPLSLTVLLL LSLLSPLQCP AA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DNASE1L1
<b>Protein Names</b>	Recommended name: Deoxyribonuclease-1-like 1 EC= 3.1.21.- Alternative name(s): DNase X Deoxyribonuclease I-like 1 Short name= DNase I-like 1 Muscle-specific DNase I-like XIB
<b>Expression Region</b>	19-302
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the deoxyribonuclease family and the protein shows high sequence similarity to lysosomal DNase I. Alternate transcriptional splice variants, encoding the same protein, have been characterized.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.