



Recombinant Schizosaccharomyces pombe Casein kinase II subunit beta (ckb1)

Product Code	CSB-YP006075SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40232
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MQLYSSESES DDSQYWVDWF LGLKGNEFFC EVDEDFIQDR FNLTGLSHEV PHYSQSLDLI LDVLDLPE EVQDEVEASA RHLYGLIHAR YILTAQGLYK MLEKYKKCDF GHCPRVLCNG QPMLPVGLSD IAHTKSVKLY CPRCEDVYTP KSQRHASIDG AYFGTSFPHM LFQVYPELAV PKSQERYIPR IFGFKVHSYS ATFKKQDVYK EKQKKRLQGA EAESKNKLAI T
Source	Yeast
Target Names	ckb1
Protein Names	Recommended name: Casein kinase II subunit beta Short name= CK II beta
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.