



Recombinant Schizosaccharomyces pombe Low molecular weight phosphotyrosine protein phosphatase (stp1)

Product Code	CSB-BP001176SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41893
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTKNIQVLFV CLGNICRSPM AEAVFRNEVE KAGLEARFDT IDSCGTGAWH VGNRPDPRTL EVLKNGIHT KHLARKLSTS DFKNFDYIFA MDSSNLRNIN RVKPQGSRAK VMLFGEYASP GVSKIVDDPY YGGSDGFGDC YIQLVDFSQN FLKSIA
Source	Baculovirus
Target Names	stp1
Protein Names	Recommended name: Low molecular weight phosphotyrosine protein phosphatase EC= 3.1.3.48 Alternative name(s): Low molecular weight cytosolic acid phosphatase EC= 3.1.3.2 PTPase Small tyrosine phosphatase
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.