



Recombinant Arabidopsis thaliana COP9 signalosome complex subunit 1 (CSN1)

Product Code	CSB-MP009858DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P45432
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MERDEEASGP MMEMCTNGGE ETSNRRPIIS GEPLDIEAYA ALYKGR TKIM RLLFIANHCG GNHALQFDAL RMAYDEIKKG ENTQLFREVV NKIGNRLGEK YGM DLAWCEA VDRRAEQKKV KLENELSSYR TNLIKESIRM GYNDFGDFYY ACGMLGDAFK NYIRTRDYCT TTKHIIHMCM NAILVSIEMG QFTHVTSYVN KAEQNPETLE PMVNAKL RCA SGLAHLELKK YKLAARKFLD VNP ELGNSYN EVIAPQDIAT YGGLCALASF DRSELKQKVI DNINFRNFLE LVPDVRELIN DFYSSRYASC LEYLASL KSN LLLDIHLHDH VDTLYDQIRK KALIQYTLPF VSVDL SRMAD AFKTSVSGLE KELEALITDN QIQARIDSHN KILYARHADQ RNATFQKVLQ MGNEFDRDVR AMLLRANLLK HEYHARSARK L
Source	Mammalian cell
Target Names	CSN1
Protein Names	Recommended name: COP9 signalosome complex subunit 1 Short name= CSN complex subunit 1 Alternative name(s): Constitutive photomorphogenesis protein 11 Protein FUSCA 6
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.