



Recombinant Meleagris gallopavo Guanine nucleotide-binding protein subunit alpha-11 (GNA11)

Product Code	CSB-EP009583MQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P45645
Product Type	Recombinant Protein
Immunogen Species	Meleagris gallopavo (Wild turkey)
Purity	>85% (SDS-PAGE)
Sequence	MTLESMMACC LSDEVKESKR INAEIEKQLR RDKRDARREL KLLLLGTGES GKSTFIKQMR IIHGSGYSEE DKKGFTKLVY QNIFTAMQSM IRAMETLKIL YKYEQNKANA VLIREVDVEK VMTFEQPYVS AIKTLWNDPG IQECYDRRRE YQLSDSAKYY LSDVDRIATP GYLPTQQDVL RVRVPTTGII EYPFDLENII FRMVDVGGQR SERRKWIHCF ENVTSIMFLV ALSEYDQVLV ESDNENRMEE SKALFRTIIT YPWFQNSSVI LFLNKKDLLE DKILYSHLVD YFPEFDGPQR DAQAAREFIL KMFVDLNPDS DKIIYSHFTC ATDTENIRFV FAAVKDTILQ LNLKEYNLV
Source	E.coli
Target Names	GNA11
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit alpha-11 Short name= G alpha-11 Short name= G-protein subunit alpha-11
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.