



Recombinant *Anas platyrhynchos* Acyl-CoA-binding protein (DBI)

Product Code	CSB-MP006519BZD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P45882
Product Type	Recombinant Protein
Immunogen Species	<i>Anas platyrhynchos</i> (Mallard) (<i>Anas boschas</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MFQAHLLRGT LTLSFFLHQA DFDEAAEEVK KLKTRPTDEE LKELYGFYKQ ATVGDINIEC PGMLDLKGKA KWEAWNLLKKG ISKEDAMNAY ISKAKTMVEK YGI
Source	Mammalian cell
Target Names	DBI
Protein Names	Recommended name: Acyl-CoA-binding protein Short name= ACBP Alternative name(s): Diazepam-binding inhibitor Short name= DBI Endozepine Short name= EP
Expression Region	1-103
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.