



Recombinant *Xenopus laevis* F-actin-capping protein subunit alpha-1 (capza1)

Product Code	CSB-EP004510XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P25229
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	EVFNDVRLLL NNDNLLREGA AHAFQAQYNMD QFTPAKIEGY DDQVLITEHG DLGNSRFLDP RNRITFKFDH LRKEASDHPH DDSDVALKSW RDACDLALRA YVKEHYPNGV CTVYGKTIDG QQTIVSCIES HQFQPKNFWN GRWRSEWKFT ISGSTAQLVG VLKIQVHYE DGNVQLVSHK DVQESITISG EAQTAKEFVK IIEQAESDYQ TAISENYQTM SDTTFKALRR QLPVTRTKID WNKILSYKIG KEMQNA
Source	E.coli
Target Names	capza1
Protein Names	Recommended name: F-actin-capping protein subunit alpha-1 Alternative name(s): Actin-binding protein chain A Short name= ABP-A
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.