



Recombinant *Saccharomyces cerevisiae* Proteasome assembly chaperone 2 (ADD66)

Product Code	CSB-MP018922SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36040
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSCLVLPLVS VGNIPQLSID WLLNSQANEW EYLEALDSKY LVEFVGPLDR PEDGSDSLYK DADMKYSSAL EVFYNNKRG LFAIQQRTPLV SVNYLNNFIV EII LPFLSKY NISEICIWDS LYAMEDENG VIVRPQEVYSL GEFYFDDEAE LLSNLHLNDQ ESMVNNWLHF TPTSFQDKIS VDQPIFKILF QILNASQRPK ALRSIKYCSC LANEGDNSLD SQQFLQWIIS QKVIKNAPPI VKFVRPISWQ GAYGMADARD KFVDLYN
Source	Mammalian cell
Target Names	ADD66
Protein Names	Recommended name: Proteasome assembly chaperone 2 Alternative name(s): Alpha-1-proteinase inhibitor-degradation deficient protein 66 Proteasome biogenesis-associated protein 2
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.