



# Recombinant Mesocricetus auratus G2/mitotic-specific cyclin-B1 (CCNB1)

<b>Product Code</b>	CSB-MP004806MRG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P37882
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mesocricetus auratus (Golden hamster)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MALRVTRNTK LNTENKAKVS MAGAKAVPVT LAAASKPGLR PRTALGDIGN KVSEQAQARL PLKKEKTSV TGKVS AKIPP PKPLEKVPPV SEPEVELAET HEPEPVMDEK LSPEPILVDN PSPSPMETSG CAPAE EYLCQ AFSDVILAVS DVDADDGADP NLCSEYVKDI YAYLRQLEEE QSVRPKYLLG REVTGNMRAI LIDWLIQVQM KFRLLQETMY MTVSIIDRFM QDNCVPPKML QLVGVTAMFI ASKYEEMYPP EIGDFAFVTN NTYTKHQIRQ MEMKILRVLN FSLGRPLPLH FLRRTSKIGE VDVEQHTLAK YLMELTLLDY DMVDFAPSQI AAGAFCLALK ILDNGEWTP T LQHLYSYTEE SLLPVMQHLA KNVVMVNHGL TKHMTIKNKY ATSKHAKIST LAQLNCTLVQ NLSKAVAKA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CCNB1
<b>Protein Names</b>	Recommended name: G2/mitotic-specific cyclin-B1
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.