



Recombinant *Saccharomyces cerevisiae* Box C/D snoRNA protein 1 (BCD1)

Product Code	CSB-MP027113SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38772
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAVLCGVCGI KEFKYKCPRC LVQTCSECS KKHKTRDNCS GQTHDPKEYI SSEALKQADD DKHERNAYVQ RDYNYLTQLK RMVHVQKMDA RMKNKRVLGP VGGHNSNFKK RRYDIDEDDR DSTECQRIIR RGVNCLMLPK GMQRSSQNRS KWDKTMDLFV WSVEWILCPM QEKGEKKELF KHVSHRIKET DFLVQQMGKN VFQKCCFYR LAGTSSCIEG EDGSETKEER TQILQKSGLK FYTKTFPYNT THIMDSKMLV ELAIHEKCIG ELLKNTTVIE FPTIFVAMTE ADLPEGYEV L HQEPRPLEHT STLNKFIDNA REEEDAEEEDS QPTEEPVQKE TQDASDSDSD SDDDYNPGLS MDFLTA
Source	Mammalian cell
Target Names	BCD1
Protein Names	Recommended name: Box C/D snoRNA protein 1
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.