



# Recombinant *Saccharomyces cerevisiae* Pre-mRNA-processing factor 17 (CDC40)

<b>Product Code</b>	CSB-EP005007SVG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P40968
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGLVDGYDTS SDSDLNFDEG KSVHEKKNNGN LHEDTSYEPS SNNIHKRKSH FTKSELKRRR KTRKGDGPWG SWSSSDDETS QASETQKEDQ DIFVHALAED NLDSEQIEVE EVSHFYGKSE KDYQGRGYLY PPNDVDVDLR EERISFRCYL PKKVIRNYPG HPEGTTALKF LPKTGHLILS GGDNDHTIKIW DFYHDYECLR DFQGHNKPIK ALRFTEDCQS FLSSSFDRSV KIWDTETGKV KTRLHLNSTP ADVESRPTNP HEFIVGLSNS KILHYDDRVS ENQGLVQTYD HHLSSILALK YFPDGSKFIS SSEDKTVRIW ENQINVPIKQ ISDTAQHSMP FLNVHPSQNY FCAQSMNRI YSFSCLKPKYK RHPKKIFKGH SSAGYGISLA FSGDGRYICS GDSKSRFTW DWNTSRLN IKIPGNKPIT QVDWHPQETS KVICS GAAGK IYVCD
<b>Source</b>	E.coli
<b>Target Names</b>	CDC40
<b>Protein Names</b>	Recommended name: Pre-mRNA-processing factor 17 Alternative name(s): Cell division control protein 40
<b>Expression Region</b>	1-455
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.