



Recombinant *Oryza sativa* subsp. japonica 60S acidic ribosomal protein P0 (Os08g0130500, LOC_Os08g03640)

Product Code	CSB-MP020336OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41095
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AIKRTKAEK KVA YDKKLCQ LLDEYTKVLI AVADNVGSNQ LQEIRKGLRG DSIVLMGKNT LIRRCIKVHA DNTGNKEFLE LMPLLVGNVG LIFTKGD LKE VREEVAKYKV GAPARVGLVA PVDVVVPPGN TGLDPSQTSF FQVLNIPTKI NKGTV EIITP VELIKKGD KV GSSESALLAK LGIRPFSYGL VITNVYDSGS VFSPEVLDLT EDDLMEKFAS GVSMVASVSL AISYPTIAAA PHMFLNGYKN VLAVAVET EY SYPHADKI KE YLKDPSKFAV AAPVAADSGA AAPSAAKEEE KKEEPEEESD GDLGMSLFD
Source	Mammalian cell
Target Names	Os08g0130500
Protein Names	Recommended name: 60S acidic ribosomal protein P0
Expression Region	2-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.