



# Recombinant Phaeosphaeria nodorum Tubulin beta chain (tubB)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP025318EUG  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P41799   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Phaeosphaeria nodorum (strain SN15 / ATCC MYA-4574 / FGSC 10173)<br>(Glume blotch fungus) (Parastagonospora nodorum)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MREIVHLQTG QCGNQIGAAF WQTISGEHGL DGSGVYNGTS DLQLERMNVY<br>FNEASGNKFV PRAVLVDLEP GTMDAVRAGP FGQLFRPDNF VFGQSGAGNN<br>WAKGHYTEGA ELVDQVLDVV RREAEGCDCL QGFQITHSLG<br>GGTGAGMGT LISKIREEFP DRMMATFSVV PSPKVS DTVV EPYNATLSIH<br>QLVENSDETF CIDNEALYDI CMRTLKLN NP SYGDLNHLVS AVMSGVTTCL<br>RFPGQLNSDL RKLAVNMV PF PRLHFFMVGF APLTSRGAHS FRAVTVPELT<br>QQMFDPKNMM AASDFRNGRY LTCSAYFRGK VSMKEVEDQM<br>RNVQKNSSY FVEWIPNNVQ TALCSVPPRG LKMSATFVGN STSIQELFKR<br>IGDQFTAMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVS EYQQYQEASI<br>SEGEEYDEE APLEAEE |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | tubB   |
| <b>Protein Names</b>     | Recommended name: Tubulin beta chain Alternative name(s): Beta-tubulin   |
| <b>Expression Region</b> | 1-447  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |