



# Recombinant Mouse Tumor necrosis factor ligand superfamily member 10 (Tnfsf10)

<b>Product Code</b>	CSB-EP023985MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P50592
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	YMYFTNEMKQLQDNYSKIGLACFSKTDEDFWDSTDGEILNRPCLQVKRQLYQ LIEEVTLR TFQDTISTVPEKQLSTPPLPRGGRPQKVAAHITGITRRSNSALIPISKDGKTLGQ KIESW ESSRKGHSFLNHVLFNRNDELVIEQEGLYIYSQTYFRFQEAEDASKMVSCKDKV RTKQLVQ YIYKYTSYPDPIVLMKSARNSCWSRDAEYGLYSIQGGLFELKKNDRIFVSVTN EHLMDL DQEASFFGAFLIN
<b>Source</b>	E.coli
<b>Target Names</b>	Tnfsf10
<b>Protein Names</b>	Recommended name: Tumor necrosis factor ligand superfamily member 10 Alternative name(s): TNF-related apoptosis-inducing ligand Short name= Protein TRAIL CD_antigen= CD253
<b>Expression Region</b>	39-291
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Target Details</b>	This protein is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This protein preferentially induces apoptosis in transformed and tumor cells, but does not appear to kill normal cells although it is expressed at a significant level in most normal tissues. This protein binds to several members of TNF receptor superfamily including TNFRSF10A/TRAILR1, TNFRSF10B/TRAILR2, TNFRSF10C/TRAILR3, TNFRSF10D/TRAILR4, and possibly also to TNFRSF11B/OPG. The activity of this protein may be modulated by binding to the decoy receptors TNFRSF10C/TRAILR3, TNFRSF10D/TRAILR4, and TNFRSF11B/OPG that cannot induce apoptosis. The binding of this protein to its receptors has been shown to trigger the activation of MAPK8/JNK, caspase 8, and caspase 3.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.