



# Recombinant *Saccharomyces cerevisiae* DNA-directed RNA polymerase III subunit RPC4 (RPC53)

<b>Product Code</b>	CSB-YP018346SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P25441
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSNKGNGRL PSLKDSSSNG GGSAPSLKF KPKAVARKSK EEREAAASKV KLEESKRGD DKKHFNKNK RVTGAGGQQR RMAKYLNNTH VISSGPLAAG NFVSEKGLR RGFIKSEGS SSVLQKLET IDNGAESSEN EAEDDDNEG ASKSKKFFNM GKEFEARNLI EDEDDGESEK SSDVMDDEE WRSKRIEQLF PVRPVRV RHE DVETVKREIQ EALSEKPTRE PTPSVKTEPV GTGLQSYLEE RERQVNEKLA DLGLEKEFQS VDGKEAAEL ELLNADHQHI LRKLKMMNNK PERFMVFQLP TRLPAPERPA VKEEKEDMET QASDPSK KKK NIKKKDTKDA LSTRELAGKV GSIRVHKSGK LSVKIGNVVM DIGKGAETTF LQDVIALSIA DDASSAELLG RVDGKIVVTP QI
<b>Source</b>	Yeast
<b>Target Names</b>	RPC53
<b>Protein Names</b>	Recommended name: DNA-directed RNA polymerase III subunit RPC4 Short name= RNA polymerase III subunit C4 Alternative name(s): C53 DNA-directed RNA polymerase III 47 kDa polypeptide
<b>Expression Region</b>	1-422
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.