



Recombinant Pisum sativum Tubulin beta-1 chain (TUBB1)

Product Code	CSB-MP025319EWE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29500
Product Type	Recombinant Protein
Immunogen Species	Pisum sativum (Garden pea)
Purity	>85% (SDS-PAGE)
Sequence	MRQILHIQGG QCGNQIGAKF WEVVCAEHGI DPTGRYTGDS DLQLERIDVY YNEASGGRFV PRAVLMDLEP GTMDSIRSGP YGQIFRPDNF VFGQSGAGNN WAKGHYTEGA ELIDSVLDVV RKEAENCDCDCL QGFQVCHSLG GGTGSGMGTL LISKIREEYP DRMMLTFSVF PSPKVSOTVV EPYNATLSVH QLVENADEVM VLDNEALYDI CFRILKLSNP SFGDLNHLIS ATMSGVTCCCL RFPGQLNSDL RKLAVNLIPF PRLHFFMVGF APLTSRGSQQ YRTLSVPELT QQMWDAKNMM CAADPRHGRY LTASAMFRGK MSTKEVDEQM MNVQNKNSSY FVEWIPNNVK STVCDIPPTG LKMASTFIGN STSIQEMFRR VSEQFTAMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVA EYQQYQDATA DEDEYGEEEG DEEEYGQHDI
Source	Mammalian cell
Target Names	TUBB1
Protein Names	Recommended name: Tubulin beta-1 chain Alternative name(s): Beta-1-tubulin
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.