



Recombinant Mouse Segment polarity protein dishevelled homolog DVL-1 (Dvl1)

Product Code	CSB-BP007284MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51141
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAETKIIYHM DEEETPYLVK LPVAPERVTL ADFKNVLSNR PVHAYKFFFK SMDQDFGVVK EEIFDDNAKL PCFNDRVSW LVLAEGAHS AGSQGTDSHT DLPPPLERTG GIGDSRPPSF HPNVASSRDG MDNETGTESM VSHRRERARR RNRDEAARTN GHPRGDRRRD LGLPPDSAST VLSSELESSS FIDSDEEDNT SRLSSSTEQS TSSRLVRKHK CRRRKQRLRQ TDRASSFSSI TDSTMSLNII TVTLNMERHH FLGISIVGQS NDRGDGGIYI GSIMKGGAVA ADGRIEPGDM LLQVNDVNF NMSNDDAVRV LREIVSQTGP ISLTVAKCWD PTPRSYFTIP RADPVRPIDP AAWLSHTAAL TGALPRYGTS PCSSAITRTS SSSLTSSVPG APQLEEAPLT VKSDMSAIVR VMQLPDSGLE IRDRMWLKIT IANAVIGADV VDWLYTHVEG FKERREARKY ASSMLKHGFL RHTVNKITFS EQCYVFGDL CSNLASLNLN SGSSGASDQD TLAPLPHPSV PWPLGQGYPY QYPGPPPCFP PAYQDPGFSC GSGSAGSQQS EGSKSSGSTR SSHRTPGREE RRATGAGGSG SESDHTVPSG SGSTGWWERP VSQLSRGSSP RSQASAVAPG LPPLHPLTKA YAVVGGPPGG PPVRELAVP PELTGSRQSF QKAMGNPCEF FVDIM</p>
Source	Baculovirus
Target Names	Dvl1
Protein Names	Recommended name: Segment polarity protein dishevelled homolog DVL-1 Short name= Dishevelled-1 Alternative name(s): DSH homolog 1
Expression Region	1-695
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	DVL1, the human homolog of the Drosophila dishevelled gene (dsh) encodes a cytoplasmic phosphoprotein that regulates cell proliferation, acting as a transducer molecule for developmental processes, including segmentation and neuroblast specification. DVL1 is a candidate gene for neuroblastomatous transformation. The Schwartz-Jampel syndrome and Charcot-Marie-Tooth disease type 2A have been mapped to the same region as DVL1. The phenotypes of these diseases may be consistent with defects which might be



expected from aberrant expression of a DVL gene during development.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.