



# Recombinant Guinea pig Tumor necrosis factor (TNF)

<b>Product Code</b>	CSB-EP023955GU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P51435
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Cavia porcellus (Guinea pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GPQREEQFSSGPPFRPLAQTTLRSASQNDNDKPVAVHVVANQQAEEELQWL SKRANALLA NGMGLSDNQLVVPDGLYLIYSQVLFKGGCPSYLLLTHTVSRLAVSYPEKVN LLSAIKS PCQKETPEGAERKPWYEPYIYLGQVFLQKGDRLSAEVNLPQYLDFAVDSGQIYF GVIAL
<b>Source</b>	E.coli
<b>Target Names</b>	TNF
<b>Protein Names</b>	Recommended name: Tumor necrosis factor Alternative name(s): Cachectin TNF-alpha Tumor necrosis factor ligand superfamily member 2 Short name= TNF-a Cleaved into the following 6 chains: 1. Tumor necrosis factor, membrane form AI
<b>Expression Region</b>	57-234
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Target Details</b>	This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.