



# Recombinant *Saccharomyces cerevisiae* UBX domain-containing protein 1 (SHP1)

<b>Product Code</b>	CSB-BP025531SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P34223
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAEIPDETIQ QFMALTNVSH NIAVQYLSEF GDLNEALNSY YASQTDDQKD RREEAHWNRQ QEKALKQEAF STNSSNKAIN TEHVGGGLCPK PGSSQGSNEY LKRKGSTSPK PTKGSSRSGS GNNSRFMSFS DMVRGQADDD DEDQPRNTFA GGETSGLEVT DPSPNSLLK DLLEKARRGG QMGAENGF RD DEDHEMGANR FTGRGFRLGS TIDAADEVVE DNTSQSQRRP EKV TREITFW KEGFQVADGP LYRYDDPANS FYLSELNQGR APLKLLDVQF GQEVEVNVYK KLDESYKAPT RKLGGFSGQG QRLGSPIDGE SSPAEV PKNE TPAAQEQPMP DNEPKQGDTS IQIRYANGKR EVLHCNSTDT VKFLYEHVTS NANTDPSRNF TLNYAFPIKP ISNDETTLKD ADLLNSVVVQ RWA
<b>Source</b>	Baculovirus
<b>Target Names</b>	SHP1
<b>Protein Names</b>	Recommended name: UBX domain-containing protein 1 Alternative name(s): Suppressor of high-copy PP1 protein
<b>Expression Region</b>	1-423
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.