



Recombinant *Notothenia coriiceps neglecta* Tubulin beta-1 chain (tubb1)

Product Code	CSB-MP025319NHX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36221
Product Type	Recombinant Protein
Immunogen Species	<i>Notothenia neglecta</i> (Yellowbelly rockcod) (<i>Notothenia coriiceps neglecta</i>)
Purity	>85% (SDS-PAGE)
Sequence	MREIVHLQAG QCGNQIGSKF WEVISDEHGI DPTGSYHGDS DLQLDRINVY YNEASGGKYV PRAVLVDLEP GTMDSVRS GP FGQIFRPDNF VFGQSGAGNN WAKGHYTEGA ELVDSVLDVV RKEAEGCDCL QGFQLTHSLG GGTGSGMGTL LISKIREEYP DRIMNTFSVV PSPKVS DTVV EPYNATLSVH QLVENTDETF CIDNEALYDI CFRTLKLTTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMVPF PRLHFFIPGF APLTSRGGQQ YRSLTVPELT QQMFDSKNMM AACDPRHGRY LTVA AIFRGR MSMKEVDEQM LNAQNKNSSY FVEWIPNNVK TAVCDIPPRG LKMAATFIGN STAIQELFKR ISEQFTAMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVS EYQQYQDATA EEEGEFEEEG EYEDGA
Source	Mammalian cell
Target Names	tubb1
Protein Names	Recommended name: Tubulin beta-1 chain Alternative name(s): Beta-1-tubulin
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.