



# Recombinant *Saccharomyces cerevisiae* Phosphoacetylglucosamine mutase (PCM1), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP017869SVG   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P38628  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | PCM1  |
| <b>Protein Names</b>     | Recommended name: Phosphoacetylglucosamine mutase Short name= PAGM<br>EC= 5.4.2.3 Alternative name(s): Acetylglucosamine phosphomutase N-<br>acetylglucosamine-phosphate mutase   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at<br>4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the<br>contents to the bottom. Please reconstitute protein in deionized sterile water to a<br>concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final<br>concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final<br>concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,<br>storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life<br>of lyophilized form is 12 months at -20°C/-80°C.  |