



Recombinant *Saccharomyces cerevisiae* Ribosomal RNA-processing protein 8 (RRP8)

Product Code	CSB-BP020531SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38961
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MALFNVEGWS IKTKTVAFDN KTNKSSKDKK KNNRKNGLT REQKLKEETE AELKEQVEDI PSEGSVAKDI PKNQEKSDQ NETSKKRKHD EEAPLMQVKE NIEKPTKKQL TPLQQKMMAK LTGSRFRWIN EQLYTISSDE ALKLIKEQPQ LFDEYHDGFR SQVQAWPENP VDFVVDQIRY RCMKPVNAPG GLPGLKDSKE IVIADMGCGE AQLALEINNF FKNYNKKAKK YLKRRHKVHS FDLKKANERI TVADIRNVPL PDESCIVVF CLALMGTNFL DFIKEAYRIL APRGELWIAE IKSRFSDGKG NEFVDALKLM GFFHKKTFFDE NKMFTREFF KPPAEIIEER RQKLERRQKF IEVETEKEEL EKKRRKIAEG KWLLKPCIYK RR
Source	Baculovirus
Target Names	RRP8
Protein Names	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.-
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.