



Recombinant Drosophila melanogaster U1 small nuclear ribonucleoprotein A (snf)

Product Code	CSB-MP022327DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43332
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MEMLPNQTIY INNLNEKIKK EELKKSLYAI FSQFGQILDI VALKTLKMRG QAFVIFKEIG SASNALRTMQ GFPFYDKPMQ IAYSKSDSDI VAKIKGTFKE RPKKVKPPKP APGTDEKKDK KKKPSSAENS NPNAQTEQPP NQILFLTNP EETNEMMLSM LFNQFPGFKE VRLVPNRHDI AFVEFTTELQ SNAAKEALQG FKITPTHAMK ITFAKK
Source	Mammalian cell
Target Names	snf
Protein Names	Recommended name: U1 small nuclear ribonucleoprotein A Short name= U1 snRNP A Short name= U1-A Short name= U1A Alternative name(s): Sex determination protein snf
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.