



Recombinant BUD13 homolog (R08D7.1)

Product Code	CSB-YP002885CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30640
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTSKADYLKK YLSPASGDIE KKKKKKNKDK NKPSGLRLIE EDAFLSVDA KTRDIGSDEE REEIEVLKQS VKKAKVVHGF KQTFAEVDAP KVIKPEPLSP DNSPPRGKRQ RHSDNSPPR PSRKRNDSDN SPPRPSRNRH DSDKDNSPPR RRRHSDNSP PRPSRKIREE SPSARNRRSP PRTRDRHDS DNSPPRNRSR RDSNSPPRR RPSSPARRRK DDDLSPPRKS RKIEEPKKIK KEEPDSDTET SGRTLEGKRS GLQSARDLKE ESKLRAKNS KMFEEMDTSV SGRFADTVYR QKQTKKKGKD SEEDQAKKER ETKKTEELKE KYKSWNKGVA QIEDRRAQLE EMARVAAEPM ARARDDAMN AHLKEVLHAA DPMANMIQKK RRDTAIDRGE LVYPSYHGHF VPNRFGIAPG YRWDGVDRSN GFEGKLAKTE NTKTANQSEY YKSIAEYE
Source	Yeast
Target Names	R08D7.1
Protein Names	Recommended name: BUD13 homolog
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.