



Recombinant *Saccharomyces cerevisiae* Ethanolamine-phosphate cytidyltransferase (ECT1)

Product Code	CSB-EP017656SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P33412
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTVNLDPDKV WIDGCFDFTH HGHAGAILQA RRTVSKENGK LFCGVHTDED IQHNKGTPVM NSSERYEHR SNRWCSEVVE AAPYVTDPNW MDKYQCQYVV HGDDITIDAN GEDCYKLVKE MGRFKVVKRT YGVSTTEIIH RILTKKSLPP THPDYYPTTQ ELSFYVAQD AVSKHCYVFQ RDLDNVLVNG GYKFDAEDCV YVDGDFDLFH MGDIDQLRKL KMDLHPDKKL IVGITSDYS STIMTKERV LSVLSCKYVD AVIIDADATS MSQYNCEKYH IGTAVLTAAG KFSEYLTKEL IVKRVESQRE VYIARNQKKG MSI
Source	E.coli
Target Names	ECT1
Protein Names	Recommended name: Ethanolamine-phosphate cytidyltransferase EC= 2.7.7.14 Alternative name(s): CTP:phosphoethanolamine cytidyltransferase Phosphorylethanolamine transferase
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.