



Recombinant General transcription factor IIH subunit 2 (T16H12.4)

Product Code	CSB-MP010008CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34567
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDDDEQKGYT WEAGYAEGLN INDVLVEDEG GSIEKSIKY VADSKRKARL TKRPERIRLG IMRHVMIVID CSRFMTSKAM PPSRFVVVMK ALQTFLDRFF EQNPQAQIGL ITCKDRKADR LTMMTGNIHV LKESLNTLTE AFCGGDFSLQ NALQLACANL KGMPGHVSRE VVLVISALST IDPGNIYSTI ETMKRMNIRC SAIGLSAEMF VCKEMAKATK GEYSVALDPD HLQLLFSKHT LPPSSAKSSE CNAIHVGFPF HELITTRSFC VCHPDTKPIS SRGFICTQCG ARHCSIPAEC PVCKLTLVAA PQLARAFRHL QPLSAFEQIE VTRGFCYACE TRLSGEGFRC GSCQLVFCLD CDTLLHESLH VCPGCN
Source	Mammalian cell
Target Names	T16H12.4
Protein Names	Recommended name: General transcription factor IIH subunit 2 Alternative name(s): TFIIH basal transcription factor complex subunit 2
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.