



Recombinant *Saccharomyces cerevisiae* Exosome complex component RRP4 (RRP4)

Product Code	CSB-MP007890SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38792
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SEVITITKR NGAFQNSSLN SYNNTGISDD ENDEEDIYMH DVNSASKSES DSQIVTPGEL VTDDPIWMRG HGTYFLDNMT YSSVAGTVSR VNRLLSVIPL KGRYAPETGD HVVGRIAEVG NKRWKVDIGG KQHAVLMLGS VNLPGGILRR KSEDELQMR SFLKEGDLN AEVQSLFQDG SASLHTRSLK YGKLRNGMFC QVPSSLIVRA KNHTHNLPGN ITVVLGVNGY IWLKTSQMD LARDTPSANN SSSIKSTGPT GAVSLNPSIT RLEEESSWQI YSDENDPSIS NNIRQAICRY ANVIKAL AFC EIGITQQRIV SAYEASMVYS NVGELIEKNV MESIGSDILT AEKMRGNGN
Source	Mammalian cell
Target Names	RRP4
Protein Names	Recommended name: Exosome complex component RRP4 Alternative name(s): Ribosomal RNA-processing protein 4
Expression Region	2-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.