



# Recombinant *Saccharomyces cerevisiae* Exosome complex component RRP4 (RRP4)

<b>Product Code</b>	CSB-YP007890SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P38792
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SEVITITKR NGAFQNSSNL SYNNTGISDD ENDEEDIYMH DVNSASKSES DSQIVTPGEL VTDDPIWMRG HGTYFLDNMT YSSVAGTVSR VNRLLSVIPL KGRYAPETGD HVVGRIAEVG NKRWKVDIGG KQHAVLMLGS VNLPGGILRR KSEDELQMR SFLKEGDLN AEVQSLFQDG SASLHTRSLK YGKLRNGMFC QVPSSLIVRA KNHTHNLPGN ITVVLGVNGY IWLKTSQMD LARDTPSANN SSSIKSTGPT GAVSLNPSIT RLEEESSWQI YSDENDPSIS NNIRQAICRY ANVIKAL AFC EIGITQQRIV SAYEASMVYS NVGELIEKNV MESIGSDILT AEKMRGNGN
<b>Source</b>	Yeast
<b>Target Names</b>	RRP4
<b>Protein Names</b>	Recommended name: Exosome complex component RRP4 Alternative name(s): Ribosomal RNA-processing protein 4
<b>Expression Region</b>	2-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.