



Recombinant *Saccharomyces cerevisiae* Securin (PDS1)

Product Code	CSB-BP019074SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40316
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MMPANEDKEN NIVYTGNESS GINFPQTPAH LLKRSHSNIL KPPVRLDQLK RDANSNNGNT LKYIQGGKEV SPTKRLHTHA QQQGRLPLAA KDNNRSKSFIFPETSNSQSKD ADLPQLQNTL SIRKNDQLRK LSQISRSRSR ANHNDLLSNS RKLQKYGSVL GYNALPKMKS LVLKDLADSG KNEESSDDDE GNEDSESKLG KKLQSALLKQ DSSDGENELN GGLGLFNEQG GLQQLIKNST KNEQKTKNDK SDKTDDYDIE IAPQRQEPLP YVPEGYSPFQ QDDIEKLKTF NSPYKLDLED EDDTPDKVDL LPLEQIDEEG EKDETECITR NQEEGAALPL LSKNFKEVAA VPTMELVYSE EGLDPEELED LVT
Source	Baculovirus
Target Names	PDS1
Protein Names	Recommended name: Securin
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.