



Recombinant *Saccharomyces cerevisiae* Low molecular weight phosphotyrosine protein phosphatase (LTP1)

Product Code	CSB-EP001176SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40347
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTIEKPKISV AFICLGNFCR SPMAEAIKFKH EVEKANLENR FNKIDSFSGTS NYHVGESPDH RTVSICKQHG VKINHKGKQI KTKHFDEYDY IIGMDESININ NLKKIQPEGS KAKVCLFGDW NTNDGTVQTI IEDPWWYGDIIQ DFEYNFKQIT YFSKQFLKKE L
Source	E.coli
Target Names	LTP1
Protein Names	Recommended name: Low molecular weight phosphotyrosine protein phosphatase EC= 3.1.3.48 Alternative name(s): Low molecular weight cytosolic acid phosphatase EC= 3.1.3.2 PTPase
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.