



Recombinant Mouse Protein Wnt-10b (Wnt10b)

Product Code	CSB-BP026130MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48614
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	NE ILGLKLPGEP PLTANTVCLT LSGLSKRQLG LCLRSPDVTA SALQGLHIAV HECQHQLRDQ RWNCSALEGG GRLPHHSAIL KRGFRESAFS FSMLAAGVMH AVATACSLGK LVSCGCGWKG SGEQDRLRAK LLQLQALSRG KTFPISQPSP VPGSVPSPPG QDTWEWGGCN HDMDFGKFS RDFLDSREAP RDIQARMRIH NNRVGRQVVT ENLKRKCKCH GTSGSCQFKT CWRAAPEFRA IGAALRERLS RAIFIDTHNR NSGAFQPRLR PRRLSGELVY FEKSPDFCER DPTLGSPGTR GRACNKTSRL LDGCGSLCCG RGHNVLRQTR VERCHCRFW CCYVLCDECK VTEWVNVCK
Source	Baculovirus
Target Names	Wnt10b
Protein Names	Recommended name: Protein Wnt-10b Alternative name(s): Protein Wnt-12
Expression Region	29-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It may be involved in breast cancer, and its protein signaling is likely a molecular switch that governs adipogenesis. This protein is 96% identical to the mouse Wnt10b protein at the amino acid level. This gene is clustered with another family member, WNT1, in the chromosome 12q13 region.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.