



# Recombinant *Saccharomyces cerevisiae* Exosome complex component RRP43 (RRP43)

<b>Product Code</b>	CSB-EP007896SVG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P25359
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AESTTLETI EIHPITFPPE VLARISPELS LQRHLSLGIR PCLRKYEEFR DVAIENNTLS RYADAGNIDT KNNILGSNVL KSGKTIVITS ITGGIIEETS AAIKDLDDFG EEELFEVTKE EDIIANYASV YPVVEVERGR VGACTDEEMT ISQKLHDSIL HSRILPKKAL KVKAGVRSAN EDGTFSVLYP DELEDDTLNE TNLKMCRKWS YVLYAKIVVL SRTGPVFDLC WNSLMYALQS VKLPRAFIDE RASDLRMTIR TRGRSATIRE TYEICDQTK SVPLMINAKN IAFASNYGIV ELDPECQLQN SDNSEEEEEVD IDMDKLNTVL IADLDTEAEE TSIHSTISIL AAPSGNYKQL TLVGGGAKIT PEMIKRSLLL SRVRADDLST RFNI
<b>Source</b>	E.coli
<b>Target Names</b>	RRP43
<b>Protein Names</b>	Recommended name: Exosome complex component RRP43 Alternative name(s): Ribosomal RNA-processing protein 43
<b>Expression Region</b>	2-394
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.