



# Recombinant *Drosophila melanogaster* 26S proteasome non-ATPase regulatory subunit 7 (Mov34)

<b>Product Code</b>	CSB-EP018912DLU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P26270
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>PSQEVSVNK VIVHPLVLLS VVDHFNRMGK IGNQKRVVGV LLGCWRSKGV  LDVSNFAVP FDEDDKDKSV WFLDHDYLEN MYGMFKKVNA  RERVVGWYHT GPKLHQNDIA INELVRRYCP NSVLVIIDAK PKDLGLPTEA  YISVEEVHDD GSPTSKTFEH VPSEIGAEAA EEVGVEHLLR DIKDTTVGSL  SQKITNQLMG LKGLNAQLRD IKQYLQRVGD SKMPINHQIV YQLQDIFNLL  PDITNDQFTG TMYVKTNDQM LVVYLASMVR SIALHNLIN NKLNRDAEE  GKSDSKEAKE KNKDSKDKDN KETKDKDGKK AEEKADKGD EGGKGSRK</p>
<b>Source</b>	E.coli
<b>Target Names</b>	Rpn8
<b>Protein Names</b>	<p>Recommended name: 26S proteasome non-ATPase regulatory subunit 7  Alternative name(s): 26S proteasome regulatory subunit RPN8 26S proteasome regulatory subunit S12 Proteasome subunit p39B Proteasome subunit p40 Protein Mov34</p>
<b>Expression Region</b>	2-338
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>