



Recombinant *Saccharomyces cerevisiae* U3 small nucleolar ribonucleoprotein protein IMP3 (IMP3)

Product Code	CSB-YP011692SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32899
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVRKCLKHHEQ KLLKKVDFLE WKQDQGHRDT QVMRTYHIQN REDYHKYNRI CGDIRRLANK LSLLPPTDPF RRKHEQLLLD KLYAMGVLTT KSKISDLENK VTVSAICRRR LPVIMHRLKM AETIQDAVKF IEQGHVRVGP NLINDPAYLV TRNMEDYVTW VDNSKIKKTL LRYRNQIDDF DFS
Source	Yeast
Target Names	IMP3
Protein Names	Recommended name: U3 small nucleolar ribonucleoprotein protein IMP3 Short name= U3 snoRNP protein IMP3 Alternative name(s): Interacting with MPP10 protein 3
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.