



# Recombinant *Synechococcus elongatus* Elongation factor Tu (tuf)

<b>Product Code</b>	CSB-BP025342FPY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P33171
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MARAKFERTK PHANIGTIGH VDHGKTTLTA AITTVLAKAG MAKARAYADI DAAPEEKARG ITINTAHVEY ETGNRHYAHV DCPGHADYVK NMITGAAQMD GAILVVSAAAD GPMPQTREHI LLAQKQGVN IVVFLNKEDM VDDAELLELV ELEVRELLSS YDFPGDDIPI VAGSALQALE AIQGGASGQK GDNPWVDKIL KLMEEVDAYI PTPEREVDRP FLMAVEDVFT ITGRGTVATG RIERGSVKVG ETIEIVGLRD TRSTTVTGVE MFQKTLDEGL AGDNVGLLLR GIQKTDIERG MVLAKPGSIT PHTKFESEVY VLKKEEGGRH TPFFPGYRPQ FYVRTTDDVTG AISDFTADDG SAAEMVIPGD RIKMTVELIN PIAIEQGMRF AIREGGRTIG AGVVSILQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	tuf
<b>Protein Names</b>	Recommended name: Elongation factor Tu Short name= EF-Tu
<b>Expression Region</b>	1-409
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.