



Recombinant Extracellular superoxide dismutase [Cu-Zn] (sod-4)

Product Code	CSB-EP022399CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34461
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	E VIRARAYIFK AEAGKIPTL IGTIDFDQSG SFLKLNQSVS GLAAGKHGFH IHEKGDGTNG CLSAGGHYNP HKLSHGAPDD SNRHIGDLGN IESPASGDTL ISVSDSLASL SGQYSIIGRS VVIHEKTDDL GRGTSDQSKT TGNAGSRLAC GTIGTVEERI LETTTASLPP VTQSQPIGSS SYYYSTFYLP IILYFLLSRI L
Source	E.coli
Target Names	sod-4
Protein Names	Recommended name: Extracellular superoxide dismutase [Cu-Zn] Short name= EC-SOD EC= 1.15.1.1
Expression Region	20-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.