



Recombinant *Saccharomyces cerevisiae* Epsin-3 (ENT3)

Product Code	CSB-EP007742SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47160
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SLEDTLANM SLYDAKKYFR KAQNVFNNT EMEGKVREAT NNEPWGASST LMDQISQGT YNFREREEILS MIFRRFTEKA GSEWRQIYKA LQLLDYLIKH GSERFIDDTR NSINLIRILE TFHYIDSQGR DQGINVRTRV KALIELLSDD NKIRAERKKA RETAKKYKGV AGGSASADGS LNSKAGFTST KVHGISVSAD FDSDNEDNED GSFSQNGYND NASRATSTPG QGKQEPEDFV DFFSSESSKP SKELIQEDEK KADEEEDDDD EFSEFQSAVP VTNPANSFNL LNTSPIEGMP ATTSSMPFYN SSTTDQGKIT PAIAEPKKVD PFSSLFSTAK ASAEAPSAPK ASQAKAAASN PVSNTTALS TDQDDDDDFG EMHGGAVQQE QNTNNNHTSS KEIDLLSF
Source	E.coli
Target Names	ENT3
Protein Names	Recommended name: Epsin-3
Expression Region	2-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.