



# Recombinant Chicken Protein Wnt-4 (WNT4)

<b>Product Code</b>	CSB-EP026137CH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P49337
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SNWLYLAK LSSVGSISEE ETCEKLGKLI QRQVQMCKRN LEVMDSVRRG AQLAIEECQY QFRNRRWNCS TLDTLPVFGK VVTQGTREAA FVYAISSAGV AFAVTRACSS GELDKCGCDR TVQGGSPQGF QWSGCSDNIA YGVAFSFSFV DVRERSKGAS SNRALMNLHN NEAGRKAILN NMRVECKCHG VSGSCEFKTC WKAMPPFRKV GNVLKEKFDG ATEVEQSEIG STKVLVPKNS QFKPHTDEDL VYLDSSPDFC DHDLKNGVLG TSGRQC�NKTS KAIDGCELMC CGRGFHTDEV EVVERCSCKF HWCCSVKCKP CHRVVEIHTC R
<b>Source</b>	E.coli
<b>Target Names</b>	WNT4
<b>Protein Names</b>	Recommended name: Protein Wnt-4
<b>Expression Region</b>	23-351
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family, and is the first signaling molecule shown to influence the sex-determination cascade. It encodes a protein which shows 98% amino acid identity to the Wnt4 protein of mouse and rat. This gene and a nuclear receptor known to antagonize the testis-determining factor play a concerted role in both the control of female development and the prevention of testes formation. This gene and another two family members, WNT2 and WNT7B, may be associated with abnormal proliferation in breast tissue. Mutations in this gene can result in Rokitansky-Kuster-Hauser syndrome and in SERKAL syndrome.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.