



Recombinant Meriones unguiculatus Interleukin-10 (IL10)

Product Code	CSB-EP011580MRA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47965
Product Type	Recombinant Protein
Immunogen Species	Meriones unguiculatus (Mongolian jird) (Mongolian gerbil)
Purity	>85% (SDS-PAGE)
Sequence	SR GEYPRNESNC THFPVSQTHM LRELRAAFSQ VKTFFQKKDQ LDNILLTDSL QQDFKGYLGC QALSEMIQFY LVEVMPQAEN HGPEIKENLN SLGEKLTLR MQLRRCHRFL PCENKSKAVE QVKNDFNKLQ EKGVYKAMNE FDIFINCIEA YMTIKMKS
Source	E.coli
Target Names	IL10
Protein Names	Recommended name: Interleukin-10 Short name= IL-10 Alternative name(s): Cytokine synthesis inhibitory factor Short name= CSIF
Expression Region	19-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.