



Recombinant *Drosophila melanogaster* 60S ribosomal protein L22 (RpL22)

Product Code	CSB-EP020182DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50887
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAPTAKTNKG DTKTAAAKPA EKKAAPAAAA AKGKVEKPKA EAAKPAAAAA KNVKKASEAA KDVKAAAAAA KPAAAKPAAA KPAAASKDAG KKAPAAAAPK KDAKAAAAPA PAKAAPAKKA ASTPAAAPPA KKAAPAKAAA PAAAAPAPAA AAPAVAKPAP KPKAKAAPAP SKVVVKNVLR GKGQKKKKVS LRFTIDCTNI AEDSIMDVAD FEKYIKARLK VNGKVNNLGN NVTFERSKLN LIVSSDVHFS KAYLKYLTKK YLKKNSLRDW IRVVANEKDS YELRYFRISS NDEDDDDAE
Source	<i>E. coli</i>
Target Names	RpL22
Protein Names	Recommended name: 60S ribosomal protein L22
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.