



Recombinant Cervus elaphus Interleukin-10 (IL10)

Product Code	CSB-MP011580DRN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51746
Product Type	Recombinant Protein
Immunogen Species	Cervus elaphus (Red deer)
Purity	>85% (SDS-PAGE)
Sequence	S RDASAPSDSS CTHFSNSLPL MLRELRTAFS RVKNFFQMKD QLDSMLLTQS LLDDFKGYLG CQALSEMIQF YLEEVMPQAE NHGPEIKEHV NSLGEKCLKL RLRLRRCHRF LPCENKSKAV EQVKS VFSKL QERGVYKAMS EFDIFINYIE TYTTMKMKN
Source	Mammalian cell
Target Names	IL10
Protein Names	Recommended name: Interleukin-10 Short name= IL-10 Alternative name(s): Cytokine synthesis inhibitory factor Short name= CSIF
Expression Region	20-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.