



Recombinant Diphthamide biosynthesis protein 1 (dph-1)

Product Code	CSB-MP007129CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49958
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MATRATQIST LVEEIANNPN LKEDLKILPS NYTFEVPKTI WKIRSTESKY VALQFPEGLI MYACVIADIL EKYTGCDTVI MGDVTYGACC VDDYTAKSMG CDLLVHYGHS CLVPIQNTDG IAMLYVFNH HINLSHLIDC VKENFQ GKRL VVVSTVQFIP SLQTLRTTFN KDDSSIRIDI PQCKPLSPGE VLGCTSPRLD ASKYDAIVYL GDGRFHLESI MIHNPEIEAF QYDPYSRKL T REFYDHDLMR KNRIGSIEIA RKCTTFGLIQ GTLGRQGNL K VVEELEAQL E RKGKKFLRVL LSEIFPEKLA MFPEVDCWVQ VACPRLSIDW GTQFPKPLLY PFELAVALDN ISVPSDHWPM DYYSNDSLGP WTNNEANRP KREKRKPHIV VRTEAS
Source	Mammalian cell
Target Names	dph-1
Protein Names	Recommended name: Diphthamide biosynthesis protein 1
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.