



# Recombinant Chicken Insulin gene enhancer protein ISL-1 (ISL1)

<b>Product Code</b>	CSB-EP011846CH-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P50211
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGDMGDPPKK KRLISLCVGC GNQIHDQYIL RVSPDLEWHA ACLKCAECNQ YLDETCTCFV RDGKTYCKRD YIRLYGIKCA KCSIGFSKND FVMRARAKVY HIECFRCVAC SRQLIPGDEF ALREDGLFCR ADHDVVERAS LGAGDPLSPL HPARPLQMAA EPISARQPAL RPHVHKQPEK TTRVRTVLNE KQLHTLRTCY AANPRPDALM KEQLVEMTGL SPRVIRVWFQ NKRCKDKKRS IMMKQLQQQQ PNDKTNIQGM TGTPMVAASP ERHDGGLQAN PVEVQSYQPP WKVLSDFALQ SDIDQPAFQQ LVNFSEGGPG SNSTGSEVAS MSSQLPDTPN SMVASPIEA
<b>Source</b>	E.coli
<b>Target Names</b>	ISL1
<b>Protein Names</b>	Recommended name: Insulin gene enhancer protein ISL-1 Short name= Islet-1
<b>Expression Region</b>	1-349
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the LIM/homeodomain family of transcription factors. The encoded protein binds to the enhancer region of the insulin gene, among others, and may play an important role in regulating insulin gene expression. The encoded protein is central to the development of pancreatic cell lineages and may also be required for motor neuron generation. Mutations in this gene have been associated with maturity-onset diabetes of the young.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.